

Double Ended Probe

This double ended probe is designed for versatile use during postmortem dissection and examination. Its dual-ended design allows for controlled probing and tissue manipulation during a variety of procedures.



SPECIFICATIONS

Dimensions

- AD130: 7" L
- AD131: 8" L
- AD132: 10" L

Materials

- Malleable stainless steel

Variations

Model No.	Description	Dimensions
AD130	Probe, Double Ended	7" L
AD131	Probe, Double Ended	8" L
AD132	Probe, Double Ended	10" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI