

# Bone Chisel, Alexander Mastoid

These stainless steel bone chisels are designed for controlled cutting and preparation of bone specimens. Their durable construction makes them suitable for scoring, cutting, scraping, cleaning, and sculpting bone across a variety of dissection applications.



## SPECIFICATIONS

### Construction

- Available in Premium or Standard grades

### Materials

- Stainless steel

### Dimensions

- Standard Grade: 6.5" L
- Premium Grade: 7" L
- Cutting Edge Width:
  - AE001 & AE101: 6 mm
  - AE002: 12 mm
  - AE003: 18 mm
  - AE102: 10 mm
  - AE103: 14 mm

- **Note:** Not intended for medical applications or patient use.

## Variations

| Model No. | Description | Edge Width | Length | Grade    |
|-----------|-------------|------------|--------|----------|
| AE001     | Bone Chisel | 6 mm       | 6.5" L | Standard |
| AE002     | Bone Chisel | 12 mm      | 6.5" L | Standard |
| AE003     | Bone Chisel | 18 mm      | 6.5" L | Standard |
| AE101     | Bone Chisel | 6 mm       | 7" L   | Premium  |
| AE102     | Bone Chisel | 10 mm      | 7" L   | Premium  |
| AE103     | Bone Chisel | 14 mm      | 7" L   | Premium  |

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI