

Bone Mallet, Lead Filled with Stainless Jacket

This bone mallet features an 8-ounce lead-filled head designed to deliver accurate and effective strikes during dissection procedures. Its balanced construction provides controlled impact and reliable force when working with bone specimens and chisels.



AE104

SPECIFICATIONS

Materials

- Premium grade stainless steel

Dimensions

- 7.5" L

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI