

Bone Mallet, Premium

This bone mallet features a 9-ounce solid stainless steel head designed for controlled and effective striking during dissection procedures. Its durable construction provides reliable impact and consistent force when working with bone specimens and chisels.



AE105

SPECIFICATIONS

Materials

- Premium grade stainless steel

Dimensions

- Overall length: 9" L
- Head length: 2" L
- Head diameter: 1"

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI