


Bone Shears, Gluck

These Gluck bone shears feature a semi-box eccentric lock and curved end for controlled cutting of bone specimens. Their specialized design provides leverage and precision for easily cutting ribs and other small bones during dissection procedures.



SPECIFICATIONS

Construction

- Available in Standard or Premium grades

Materials

- Stainless steel

Dimensions

- AF001 & AF101: 8.5" L

Note:

- Not intended for medical applications or patient use.

Variations

Model No.	Description	Grade	Dimensions
AF001	Bone Shears, Gluck	Standard	8.5" L
AF101	Bone Shears, Gluck	Premium	8.5" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI