

# Rib Cutting Knife

This rib cutter is uniquely designed for controlled and efficient cutting of ribs during postmortem examinations. Its specialized construction provides reliable leverage and precision when working with bone during detailed autopsy procedures.



## SPECIFICATIONS

### Construction

- Cutting Edge: Extra thick

### Materials

- Standard grade stainless steel

### Dimensions

- Overall length:  
AF007: 7.75" L  
AF008: 8" L  
AF008-M: 9" L
- Cutting edge length:  
AF007: 1.6" L  
AF008: 0.98" L  
AF008-M: 1" L

### Note:

- Not intended for medical applications or patient use. Product may vary.

### Variations

| Model No. | Description       | Cutting Edge | Dimensions |
|-----------|-------------------|--------------|------------|
| AF007     | Rib Cutting Knife | 1.6" L       | 7.75" L    |
| AF008     | Rib Cutting Knife | 0.98" L      | 8" L       |
| AF008-M   | Rib Cutting Knife | 1" L         | 9" L       |

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI