

Bone Chisel, U.S. Army Pattern

This U.S. Army pattern bone chisel is designed for controlled cutting and preparation of bone specimens. Also known as an osteotome, its durable design allows for scoring, cutting, scraping, cleaning, and sculpting bone during a variety of dissection procedures.



SPECIFICATIONS

Construction

- U.S. Army pattern
- Available in Standard or Premium grades

Materials

- Stainless steel

Dimensions

- 6.5" L
- Cutting Edge Width:
AG002 & AG102: 6 mm
AG002-H: 6 mm Heavy Duty
AG003 & AG103: 12 mm
AG004 & AG104: 18 mm

- **Note:** Not intended for medical applications or patient use.

Variations

Model No.	Description	Edge Width	Length	Grade
AG002	Bone Chisel	6mm	6.5" L	Standard
AG002-H	Bone Chisel	6 mm Heavy Duty	6.5" L	Standard
AG003	Bone Chisel	12 mm	6.5" L	Standard
AG004	Bone Chisel	18 mm	6.5" L	Standard
AG102	Bone Chisel	6mm	6.5" L	Premium
AG103	Bone Chisel	12 mm	6.5" L	Premium
AG104	Bone Chisel	18 mm	6.5" L	Premium

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI