

Bone Chisel, Lambotte

This Lambotte bone chisel is designed for controlled cutting and preparation of bone specimens. Also known as an osteotome, its specialized design allows for scoring, cutting, scraping, cleaning, and sculpting bone during a variety of procedures.



SPECIFICATIONS

Construction

- Lambotte pattern
- Available in Standard or Premium grades

Materials

- Stainless steel

Dimensions

- Standard Grade: 8" L
- Premium Grade: 9" L
- Cutting Edge Width:
AG005 & AG105: 0.5"
AG006 & AG106: 1"

- **Note:** Not intended for medical applications or patient use.

Variations

Model No.	Description	Edge Width	Length	Grade
AG005	Bone Chisel	0.5" L	8" L	Standard
AG006	Bone Chisel	1" L	8" L	Standard
AG105	Bone Chisel	0.5" L	9" L	Premium
AG106	Bone Chisel	1" L	9" L	Premium

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI