

Bone Chisel, Brunetti's

These Brunetti chisels are specifically designed for controlled vertebral dissection and access around the spinal cord. Their specialized design allows for precise cutting and preparation of vertebral bone while providing greater manual control during detailed dissection procedures.



SPECIFICATIONS

Construction

- Right or left handed options available

Materials

- Standard grade stainless steel

Dimensions

- AG015 & AG016: 11" L

- **Note:** Not intended for medical applications or patient use.

Variations

Model No.	Description	Orientation	Dimensions
AG015	Bone Chisel	Right-handed	11" L
AG016	Bone Chisel	Left-handed	11" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI