


Cleveland Bone Rongeurs No. 4

These Cleveland Rongeurs are designed for controlled cutting, sculpting, and trimming of both soft and hard bone. Their durable construction provides reliable performance and lasting use when working with bone specimens during dissection procedures.



AG033

SPECIFICATIONS

Materials

- German grade satin stainless steel

Dimensions

- 5.5" (139 mm) L
- Jaw size: 6 mm

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI