


Bone Rongeurs, Cleveland No. 4A

These Cleveland Rongeurs are designed for controlled cutting, sculpting, and trimming of both soft and hard bone. Their satin surface helps reduce glare for improved visibility, while the durable construction provides reliable performance during dissection procedures.



AG034

SPECIFICATIONS

Materials

- German grade stainless steel

Dimensions

- 6.5" (165 mm) L
- Jaw size: 6 mm

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI