

Bone Rongeurs, Beyer

This curved, double-action bone rongeur is designed for precise removal of small bone fragments during research and dissection procedures. Its curved profile and double-action mechanism provide controlled access and efficient cutting when working with bone specimens.



AG035

SPECIFICATIONS

Materials

- German grade stainless steel

Dimensions

- 7" (177.8 mm) L
- Jaw size: 3 mm

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI