

Littauer-Liston Bone Cutting Forceps, 6"

These Littauer-Liston bone-cutting forceps are designed for controlled bone cutting during autopsy and post-mortem dissection procedures. The sharp cutting tips support clean cutting, while the screw-lock mechanism helps securely hold bone during use.



AG110

SPECIFICATIONS

Construction

- Littauer-Liston pattern
- Screw-lock mechanism
- Smooth grip
- Sharp cutting tips

Performance

- Designed for controlled bone cutting
- Screw lock helps securely hold bone
- Sharp tips support clean, precise cutting

Materials

- Premium grade stainless steel

Dimensions

- 6" L (152 mm)

Note:

- **Not intended for medical applications or patient use.**

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI