

Virchow Skull Breaker

This Virchow skull breaker features a T-handle grip and thin access blade designed for controlled separation of cranial plates during autopsy procedures. Its specialized construction provides secure handling and precise leverage to assist in separating and exposing cranial structures during postmortem examinations.



SPECIFICATIONS

Construction

- Available in Standard or Premium grades

Materials

- Stainless steel

Dimensions

- 3.625" (92 mm) L
- Cutting blade width:
AH004: 0.625" W
AH104: 0.5" W

Note:

- Not intended for medical applications or patient use.

Variations

Model No.	Description	Grade	Cutting Blade Width
AH004	Virchow Skull Breaker	Standard	0.625" W
AH104	Virchow Skull Breaker	Premium	0.5" W

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI