


Bone Mallet, Postmortem, Hook Handle

This hook handled mallet is designed for controlled striking during autopsy and postmortem dissection procedures. Its specialized handle provides a secure grip and improved control for delivering consistent, effective impact when working with bone specimens.



SPECIFICATIONS

Construction

- Available in Standard or Premium grades

Materials

- 12 oz. solid stainless steel

Dimensions

- Standard Grade: 8" L
- Premium Grade: 9.5" L

Note:

- Not intended for medical applications or patient use.

Variations

Model No.	Description	Grade	Dimensions
AH005	Bone Mallet, Hook Handle	Standard	8" L
AH105	Bone Mallet, Hook Handle	Premium	9.5" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI