

Knife, Heavy Dissecting, Optional Ruler

This extra-thick, rigid stainless steel blade is designed for controlled cutting through hard tissue. Its 0.062" (1.6 mm) thick top edge tapers into a sharp cutting edge, with an optional imprinted metric ruler for convenient measurement.



SPECIFICATIONS

Construction

- Cutting Edge: Thick tapered edge with optional 40 mm ruler (AH019)

Materials

- Standard grade stainless steel

Dimensions

- Overall length: 13" (330 mm) L
- Blade length: 8" (203 mm) L X 0.75" (20 mm) W

Note:

- Not intended for medical applications or patient use.

Variations

Model No.	Description	Type	Dimensions
AH018	Rib Cutting Knife	Without Ruler	13" L
AH019	Rib Cutting Knife	With Ruler	13" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI