


Tenaculum, Dissecting

This dissecting tenaculum features a sharp, rounded tip and secure-grip handles for controlled tissue grasping and holding during dissection procedures. Its specialized design provides a firm grip and precise handling, allowing for reliable tissue manipulation during detailed examinations.



AH031

SPECIFICATIONS

Materials

- Standard grade stainless steel

Dimensions

- 6.5" (165 mm) L

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI