

Knife, German Dissecting

This lightweight, well-balanced postmortem dissection knife features a hollow-ground cutting edge for precise, controlled cutting. Its smooth handle provides added comfort, while the hardened blade is designed to withstand repeated sharpening for lasting performance.



AH101

SPECIFICATIONS

Materials

- Premium grade stainless steel

Dimensions

- 11.75" (299 mm) L X 0.79" (20 mm) W

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI