

Knife, Large Grossing

This large grossing knife features a non-slip handle and well-balanced, honed blade for comfortable, controlled cutting. Its durable blade can be repeatedly sharpened for long-lasting performance and easily sanitized with disinfectant.



AH110

SPECIFICATIONS

Construction

- Cutting Edge: Honed with air-gap grooves

Materials

- Premium grade stainless steel with plastic handle

Dimensions

- Overall length: 19.75" (501.7 mm) L
- Blade length: 14" (355 mm) L X 1.5" (38 mm) W

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI