

Knife, Autopsy Grossing, Low Profile

This autopsy grossing knife features a non-slip handle and securely bonded, well-balanced blade for comfortable, controlled cutting. Its hand-ground, honed edge can be repeatedly sharpened for long-lasting performance and easily sanitized with disinfectant.



SPECIFICATIONS

Materials

- Premium grade stainless steel

Dimensions

- Blade length:
AH114: 6" L
AH115: 7" L

Note:

- Not intended for medical applications or patient use.

Variations

Model No.	Description	Dimensions
AH114	Knife, Autopsy Grossing	6" L
AH115	Knife, Autopsy Grossing	7" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI