

Scalpel Blades, Double-Edged, 100 pack

These double-edged scalpel blades are designed for trimming and sectioning specimens during post-mortem and laboratory procedures. The dual cutting edges provide a versatile cutting surface for routine specimen preparation.



AJ111

SPECIFICATIONS

Construction

- Double-edged blade
- Non-sterile
- Slotted center opening

Material

- Stainless steel

Performance

- Designed for specimen trimming
- Dual cutting edges support versatile laboratory use

Quantity

- 100 per pack

Dimensions

- 2.25" X 0.75" (57 X 19 mm)
- Inside slot: 1.25" (32 mm)

Note:

- Not intended for medical applications or patient use.

MOPEC
GROUP

Mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI