

# Stainless Steel #3L Scalpel Handle

This #3L scalpel handle is designed for secure, controlled handling during post-mortem and laboratory procedures. The extended handle design provides additional reach, while the textured grip areas support a stable grasp during use.



AJ142

## SPECIFICATIONS

### Performance

- Designed for secure, controlled handling
- Extended design provides additional reach
- Fits #10 through #15 scalpel blades

### Construction

- #3L handle
- Extended design
- Textured grip areas

### Material

- Stainless steel

### Note:

- Not intended for medical applications or patient use.

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI