

Standard Grade #20 Scalpel Blades, 100 Pack

These #20 scalpel blades are designed for controlled cutting during post-mortem and laboratory procedures. The broad, curved cutting edge supports smooth, precise sectioning, while the standard-grade carbon steel construction provides dependable performance for routine use.



AJ143

SPECIFICATIONS

Construction

- #20 blade
- Curved cutting edge
- Non-sterile
- Individually wrapped

Material

- Standard grade carbon steel

Performance

- Designed for controlled cutting
- Broad blade profile supports smooth sectioning
- Fits #4 #4L, and #6 scalpel handles

Quantity

- 100 per pack

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI