

#10 Scalpel Blades, Sterile, 100 Pack

These sterile #10 scalpel blades are designed for controlled cutting during post-mortem and laboratory procedures. The curved cutting edge supports smooth incisions and soft-tissue cutting, while the traditional blade profile provides versatile use.



AJ170

SPECIFICATIONS

Construction

- #10 blade
- Curved cutting edge
- Sterile
- Individually wrapped

Material

- Carbon steel

Performance

- Designed for controlled cutting
- Suitable for incisions and soft-tissue cutting
- Fits #3, #3L, and #5 scalpel handles

Quantity

- 100 per pack

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI