

#12B Scalpel Blades, Sterile, 100 Pack

These sterile #12B scalpel blades are designed for precise cutting during post-mortem and laboratory procedures. The small curved blade is sharpened on both sides of the curved edge, supporting controlled specialty cutting.



AJ173

SPECIFICATIONS

Construction

- #12B blade
- Small curved blade
- Both sides of curved edge sharpened
- Sterile
- Individually wrapped

Material

- Carbon steel

Performance

- Designed for precise, controlled cutting
- Suitable for specialty cutting procedures
- Fits #3, #3L, and #5 scalpel handles

Quantity

- 100 per pack

MOPEC
GROUP

Mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI