

#15C Scalpel Blades, Sterile, 100 Pack

These sterile #15C scalpel blades are designed for precise cutting during post-mortem and laboratory procedures. The narrow curved blade and sharp point support controlled, detailed cutting, while the longer cutting edge provides additional reach compared with a #15 blade.



AJ175

SPECIFICATIONS

Construction

- #15C blade
- Narrow curved blade
- Sharp point
- Longer cutting edge
- Sterile
- Individually wrapped

Material

- Carbon steel

Performance

- Designed for precise, controlled cutting
- Suitable for detailed specialty procedures
- Fits #3, #3L, and #5 scalpel handles

Quantity

- 100 per pack

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI