

Stainless Steel #11 Scalpel Blades, 100 Pack

These sterile #11 scalpel blades are designed for precise cutting during post-mortem and laboratory procedures. The triangular blade and sharp pointed tip support small, controlled incisions and detailed cutting.



AJ204

SPECIFICATIONS

Construction

- #11 blade
- Triangular blade profile
- Sharp pointed tip
- Sterile
- Individually wrapped

Material

- Stainless steel

Performance

- Designed for precise, controlled cutting
- Suitable for small or stab-style incisions
- Fits #3, #3L, and #5 scalpel handles

Quantity

- 100 per pack

MOPEC
GROUP

Mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI