

Stainless Steel #15 Scalpel Blades, 100 Pack

These sterile #15 scalpel blades are designed for precise cutting during post-mortem and laboratory procedures. The small, curved cutting edge supports short, controlled incisions and detailed cutting.



AJ207

SPECIFICATIONS

Construction

- #15 blade
- Small curved cutting edge
- Sterile
- Individually wrapped

Material

- Stainless steel

Performance

- Designed for short, precise incisions
- Suitable for detailed, controlled cutting
- Fits #3, #3L, and #5 scalpel handles

Quantity

- 100 per pack

MOPEC
GROUP

Mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI