

Stainless Steel #20 Scalpel Blades, 100 Pack

These sterile #20 scalpel blades are designed for controlled cutting during post-mortem and laboratory procedures. The broad, curved cutting edge supports smooth, precise sectioning and routine specimen preparation.



AJ209

SPECIFICATIONS

Construction

- #20 blade
- Broad curved cutting edge
- Sterile
- Individually wrapped

Material

- Stainless steel

Performance

- Designed for controlled cutting
- Supports smooth sectioning and specimen preparation
- Fits #4, #4L, and #6 scalpel handles

Quantity

- 100 per pack

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI