

Ruler, Foldable

This large, foldable 8-foot ruler features inch and centimeter markings for accurate measurements across a variety of applications. Its convenient design folds in half for compact storage and includes 0.25" diameter holes at both ends for easy wall hanging when not in use.



BA048

SPECIFICATIONS

Materials

- Stainless steel with plastic faced ruled edge

Dimensions

- 8' (243 cm) L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI