

Evidence Bag Heat Sealer

A heat sealing device designed for securely closing evidence bags and maintaining chain of custody. It creates a strong, tamper-evident seal on plastic evidence packaging, helping ensure integrity during storage and transport in forensic and law enforcement applications.



SPECIFICATIONS

Construction

- Temperature control knob

Materials

- Heavy-duty Kraft paper stock

Performance

- Seals bags in the field or in the lab
- No warm-up time
- Uses 120V AC current

Variations

Model No.	Description	Size
FS158	Evidence Bag Heat Sealer	8"
FS159	Evidence Bag Heat Sealer	12"
FS160	Evidence Bag Heat Sealer	20"

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI