


Multi-Barrier Pouch

The Multi-barrier Pouch is a tamper-evident, seven-layer laminated pouch designed to protect FTA® cards from exposure to gases and liquids. It includes an outer paper layer for labeling and comes in small and large sizes.



SPECIFICATIONS

Construction

- Tamper-evident seal

Quantity

- 100 per package

Performance

- Outer paper layer can be used for labeling or writing
- Maintains sample security
- Protect the card from exposure to any gaseous or liquid external element

Materials

- (7) laminated layers

Variations

Model No.	Description	Dimensions
FS161	Multi-Barrier Pouch	3.7" L X 3" W
FS162	Multi-Barrier Pouch	4.4" L X 6.5" W

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI