

Probe, Flexible with Eyelet

This flexible probe features an eyelet for thread or string tracking during postmortem procedures. Its flexible design allows for controlled navigation through tissue and anatomical structures.



SPECIFICATIONS

Dimensions

- AD014: 6" L
- AD015: 10" L
- AD026: 5.5" L
- AD027: 8" L
- AD028: 12" L
- AD133: 7" L

Materials

- Standard grade stainless steel

Quantity

- 6 per pack

- **Note:** Not intended for medical applications or patient use

Variations

Model No.	Description	Dimensions
AD014	Probe, Flexible with Eyelet	6" L
AD015	Probe, Flexible with Eyelet	10" L
AD026	Probe, Flexible with Eyelet	5.5" L
AD027	Probe, Flexible with Eyelet	8" L
AD028	Probe, Flexible with Eyelet	12" L
AD133	Probe, Flexible with Eyelet	7" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI