

# Rochester Ochsner Hemostatic Forceps

Rochester Ochsner hemostatic forceps are designed for secure tissue clamping during post-mortem procedures. Serrated 1 x 2 toothed tips provide a firm, controlled grip, with straight or curved configurations available in standard or premium options.



## SPECIFICATIONS

### Construction

- Tips: Serrated, 1 x 2 teeth
- Grip: Finger ring
- Adjustable interlocking cross arms
- Available in straight or curved configurations

### Performance

- Provides secure, controlled clamping
- 1 x 2 teeth provide additional tissue grip
- Curved options provide additional access for varied laboratory applications

### Dimensions

- Available in 6.25", 7.25", 8", 9", or 10" L

### Materials

- Standard grade stainless steel

### Note:

- **Not intended for medical applications or patient use.**

## Variations

| Model No. | Description                | Edge     | Length |
|-----------|----------------------------|----------|--------|
| AC017     | Forceps, Rochester Ochsner | Straight | 6.25"  |
| AC018     | Forceps, Rochester Ochsner | Straight | 7.25"  |
| AC019     | Forceps, Rochester Ochsner | Straight | 8"     |
| AC020     | Forceps, Rochester Ochsner | Straight | 9"     |
| AC021     | Forceps, Rochester Ochsner | Straight | 10"    |
| AC022     | Forceps, Rochester Ochsner | Curved   | 6.25"  |
| AC023     | Forceps, Rochester Ochsner | Curved   | 7.25"  |
| AC024     | Forceps, Rochester Ochsner | Curved   | 8"     |
| AC025     | Forceps, Rochester Ochsner | Curved   | 9"     |
| AC026     | Forceps, Rochester Ochsner | Curved   | 10"    |

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI