

CUTTING BOARD IS TO
BE FREE FROM SCRATCHES,
STAINS, DISCOLORATION, &
DEFORMATION

REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM







MOPEC LEGAL NOTE
THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC MANAGER
MOPEC 2025

UNLESS OTHERWISE SPECIFIED:
ALL DIMENSIONS IN INCHES
TOLERANCES
0.0 ± 0.1 in
0.00 ± 0.05 in
0.000 ± 0.010 in
0.0000 ± 0.0010 in
ANGLES ±0.5°

BREAK ALL CORNERS 0.063 X 45°
ALL FILLET RADII R0.063 /0.125in
INSPECT OVAL DIMENSIONS

MATERIAL:


MOPEC GROUP
Mopec MORTECH MANUFACTURING
DRAWN TO 3RD ANGLE PROJECTION

DRAWN BY: jmarshall
DATE: 7/1/2025
CHECKED BY: ---
DATE: 7/1/2025
APPROVED BY: M.McCLAIN
DATE: 7/1/2025

PART NUMBER: AX065-096-A
DESCRIPTION: DISSECTION BOARD, WHITE W/ 3" LEGS, 96" GROSSING STATION
SIZE: A2
FILE NAME: AX065-096-A.dft
REV: A
SCALE: 1:2
DO NOT SCALE IF IN DOUBT ASK
SHEET 1 OF 3

AX065-096-A.dft

REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM



Item #	Qty.	File Name	Title	Material	Rev.
1	8	PB0609.par	TAPPING INSERT, 1/4"-20 INTERNAL, BRASS	Stainless steel	A
2	8	PB0611.par	SCREW, 1/4"-20 x 3/4"L, PAN, 18-8 SS	McMASTER-CARR # 91772A540	A
3	8	PM0618.par	Rubber Bumper, 13/16" OD, 13/16" High		A
4	8	S6829-B.asm	3" LEG ASSEMBLY		B
5	2	S7540-A.par	DISSECTION BOARD, 21" x 18.5", WHITE	3/4" HDPE - NATURAL, TEXTURED	A

MOPEC LEGAL NOTE
THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC. MANAGER
MOPEC 2025

UNLESS OTHERWISE SPECIFIED:
ALL DIMENSIONS IN INCHES
TOLERANCES
0.0 ± 0.1 in
0.00 ± 0.05 in
0.000 ± 0.010 in
0.0000 ± 0.0010 in
ANGLES ±0.5°
BREAK ALL CORNERS 0.063 X 45°
ALL FILLET RADII R0.063 /0.125in
INSPECT OVAL DIMENSIONS

MATERIAL:

MOPEC GROUP
Mopec MORTECH MANUFACTURING

DRAWN TO 3RD ANGLE PROJECTION


DRAWN BY: lmarshall 7/11/2025
CHECKED BY: --- 7/11/2025
APPROVED BY: M.McCLAIN 7/11/2025

PART NUMBER: AX065-096-A
DESCRIPTION: DISSECTION BOARD, WHITE W/ 3" LEGS, 96" GROSSING STATION
SIZE: A2 FILE NAME: AX065-096-A.dft REV: A
SCALE: 1:2 DO NOT SCALE IF IN DOUBT ASK
SHEET 2 OF 3


AX065-096-A.dft


CUTTING BOARD IS TO
BE FREE FROM SCRATCHES,
STAINS, DISCOLORATION, &
DEFORMATION


REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM



BOTTOM







MOPEC LEGAL NOTE

THIS DOCUMENT CONTAINS
PROPRIETARY INFORMATION OF
MOPEC, AND MAY NOT BE COPIED IN WHOLE
OR PART, AND NO USE SHALL BE MADE OF
IT OTHER THAN FOR WHICH IT HAS BEEN
SUBMITTED WITHOUT THE WRITTEN
PERMISSION OF A MOPEC MANAGER

MOPEC 2025

UNLESS OTHERWISE SPECIFIED
ALL DIMENSIONS IN INCHES
TOLERANCES
0.0 ± 0.1 in
0.00 ± 0.05 in
0.000 ± 0.010 in
0.0000 ± 0.0010 in
ANGLES ±0.5°

BREAK ALL CORNERS 0.063 X 45°
ALL FILLET RADII R0.063 /0.125in
INSPECT OVAL DIMENSIONS

MATERIAL:

MOPEC GROUP
Mopec MORTECH MANUFACTURING

DRAWN TO
3RD ANGLE
PROJECTION

DRAWN BY: jmarshall
DATE: 7/11/2025
CHECKED BY: ---
DATE: 7/11/2025
APPROVED BY: MMcCLAIN
DATE: 7/11/2025

PART NUMBER: AX065-096-A
DESCRIPTION: DISSECTION BOARD, WHITE W/ 3" LEGS, 96" GROSSING STATION
SIZE: A2
FILE NAME: AX065-096-A.dft
REV: A
SCALE: 1:2.67
DO NOT SCALE IF IN DOUBT ASK
SHEET 3 OF 3

AX065-096-A.dft