KITS

Student Instrument Case

This student instrument case is designed to securely store and organize up to twelve basic post-mortem dissection instruments. The zippered vinyl exterior and velvet-lined interior provide a compact, discreet storage solution for laboratory and instructional use.



AA300



AA301

SPECIFICATIONS

Construction
- Zippered vinyl case
- Dark red velvet-lined interior

Materials
- Vinyl

Performance
- Holds up to 12 basic post-mortem dissection instruments
- Supports secure, organized instrument storage

Dimensions
- 7" W × 8.25" H (178 × 209.5 mm), closed
- 15" (381 mm) end-to-end, open

Note:
- **Not intended for medical applications or patient use.**

Variations

Model No.	Description	Color
AA300	Instrument Case	Black
AA301	Instrument Case	Red

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI