BLADE REMOVER

Aunex Scalpel Blade Remover

The Aunex ABR-100U provides single-handed scalpel blade removal using a simple push-and-pull motion. The system securely disengages and contains used blades without requiring fingers or forceps and accommodates both standard and larger pathology blade sizes.



ABR-100U

SPECIFICATIONS

Performance

- Supports single-handed blade removal
- Removes blades without fingers or forceps
- Push-and-pull operation
- Fits standard blade sizes #6 through #26
- Fits larger #40, #60, and #70 blades
- Holds up to 200 blades (based on #10 blades)

Construction

- Uniflex blade-removal technology
- Lifts and traps the heel of the blade before disengaging it
- Enclosed blade containment system

MOPEC GROUP

Mopec

MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI