RETRACTORS AND LIFTERS

Double-Ended Section Lifter

This double-ended section lifter is designed to separate and lift tissue during post-mortem dissection procedures. Its stainless steel construction provides a durable instrument for routine dissection use.



AD040

SPECIFICATIONS

Construction
- Double-ended section lifter
- Flat lifting ends

Materials
- Standard-grade stainless steel

Performance
- Designed to separate and lift tissue
- Double-ended design supports versatile use during dissection

Dimensions
- Overall length: 6”

Note:
- **Not intended for medical applications or patient use.**

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI