RETRACTORS AND LIFTERS

Light-Duty Dissection Tool

This light-duty dissection tool is designed for delicate post-mortem procedures, including tissue isolation, retraction, and teasing. Its narrow dual-prong tip supports controlled handling, while the serrated grip provides a secure hold during use.



AD120

SPECIFICATIONS

Construction

- Dual-prong
- Sharp, narrow tip
- Serrated perimeter grip

Performance

- Designed for delicate, light-duty dissection
- Supports tissue isolation, retraction, and teasing

Dimensions

- Overall length: 6.125” (169 mm)

Note:

- **Not intended for medical applications or patient use.**

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI