RETRACTORS AND LIFTERS

Dual-Prong Dissection Tool

This dual-prong dissection tool is designed for delicate, light-duty post-mortem procedures. Its two sharp teeth support controlled tissue isolation, retraction, and teasing, while the serrated grip helps provide secure handling during use.



AD122

SPECIFICATIONS

Construction

- Dual-prong
- Two sharp 5 mm teeth
- Serrated perimeter grip

Performance

- Designed for delicate, light-duty dissection
- Supports tissue isolation, retraction, and teasing

Dimensions

- Overall length: 6.25” (159 mm)

Note:

- **Not intended for medical applications or patient use.**

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI