SCALPEL HANDLES

Premium Ergonomic Scalpel Handles

These premium ergonomic scalpel handles are designed for comfortable handling and are available in three sizes to accommodate different scalpel blade ranges commonly used in anatomical pathology.



AJ108



AJ109



AJ110

SPECIFICATIONS

Performance

- Ergonomic design for comfortable handling
- Available in three handle sizes
- Fits a range of commonly used scalpel blades

Compatibility

- #5 Handle fits #6 – #16 blades
- #6 Handle fits #20 – #26 blades
- #8 Handle fits #60 – #70 blades

Note:

- **Not intended for medical applications or patient use.**

Variations

Model No.	Description	Handle Size	Blade Range
AJ108	Disposable Scalpels	#5	#6-#16 blades
AJ109	Disposable Scalpels	#6	#20-#26 blades
AJ110	Disposable Scalpels	#8	#60-#70 blades

MOPEC GROUP

Mopec

MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI