SCALPEL BLADES

Stainless Steel #10 Scalpel Blades, 100 Pack

These sterile #10 scalpel blades are designed for controlled cutting during post-mortem and laboratory procedures. The curved cutting edge supports smooth incisions and soft-tissue cutting, while the traditional blade profile provides versatile use.



AJ203

SPECIFICATIONS

Construction
- #10 blade
- Curved cutting edge
- Sterile
- Individually wrapped

Material
- Stainless steel

Performance
- Designed for controlled cutting
- Suitable for incisions and soft-tissue cutting
- Fits #3, #3L, and #5 scalpel handles

Quantity
- 100 per pack

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI