
GROSSING STATION ACCESSORIES

Ensemble CPU Mount



SPECIFICATIONS

Construction
- Designed to attach to the top of your Mopec Maestro grossing station
- Width adjustable by loosening/ tightening wingnut by hand
- Phillips Screwdriver required for installation

Capacity
- 20 lbs

Performance
- Can go on either the left or right side of the station
- Mounts high on the station to be out of the user's way, and out of the way of pegboard wings and other pegboard accessories
- Accommodates any CPU that is 2.5" - 5.25" wide

Materials
- 304 stainless steel construction with a polished No. 4 finish

Dimensions
- 12" L x 8" D x 12.5" H

Compatibility
- Fits all Mopec Maestro and Maestro Encore grossing stations

MOPEC GROUP

Mopec ELEVATING PATHOLOGY

MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI