Instrument Jar, Weighted Base

This instrument holder is designed to securely accommodate a variety of pathology tools, including forceps, scissors, scalpels, and rulers. Its heavily weighted base helps prevent tipping during use, while the easy-to-clean, autoclavable design allows for convenient maintenance and sterilization.



SPECIFICATIONS

Dimensions

- BA008: 4.5" H
- BA009: 7.25" H
- BA111: 4" H
- Diameter:
  BA008 & BA009: 2.25"
  BA111: 1"

Materials

- Stainless steel

Performance

- Autoclavable

Variations

Model No.	Description	Diameter	Dimensions
BA008	Weighted Base Jar	2.25"	4.5" H
BA009	Weighted Base Jar	2.25"	7.25" H
BA111	Weighted Base Jar	1"	4" H

MOPEC GROUP
Mopec ELEVATING PATHOLOGY
MORTECH MANUFACTURING
www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI