GLOVES AND GLOVE DISPENSERS

Stainless Steel Glove Dispenser

This adjustable glove dispenser is designed to hold up to three glove boxes while accommodating different box sizes. The upper section slides vertically for height adjustment and the dispenser can be mounted using double-sided foam tape or hardware.



BA027

SPECIFICATIONS

Construction

- Sliding upper section for height adjustment
- Mounting holes for hardware installation
- Can also mount with double-sided foam tape

Performance

- Holds up to 3 glove boxes
- Adjustable height accommodates up to three boxes
- Designed for wall mounting
- Gloves not included

Dimensions

- 5" W × 4" D
- Adjustable height: 8.5"–15"

Materials

- 16-gauge cold-rolled steel

MOPEC GROUP
Mopec
MORTECH MANUFACTURING
www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI