Instrument Sterilization Trays

This three-piece stainless steel tray set is designed for sterilization and decontamination of pathology dissection instruments during laboratory procedures. Available in multiple sizes and volume capacities, it provides a convenient solution for instrument cleaning.



BA054-2H



BA054

SPECIFICATIONS

Construction
- Three pieces include instrument tray, perforated insert tray, and lid

Materials
- 18/8, 22 gauge stainless steel
- 18% chromium to 8% nickel

Options & Accessories
- For use with Instrument Cleaner (BE046)

Dimensions
- BA054: 20.75" L X 12.75" W X 6" D
- BA054-2H: 12.75" L X 10.375" W X 2" D
- BA054-4: 20.75" L X 12.75" W X 4" D
- BA054-4H: 12.75" L X 10.375" W X 4" D
- BA054-6H: 12.75" L X 10.375" W X 6" D

Capacity
- BA054: 22 qt
- BA054-2H: 4 qt
- BA054-4: 15 qt
- BA054-4H: 7 qt
- BA054-6H: 10 qt

Performance
- Autoclavable
- Resistant to corrosion and oxidation

Variations

Model No.	Description	Capacity	Dimensions
BA054	Instrument Sterilization Trays	22 qt	20.75" L X 12.75" W X 6" D
BA054-2H	Instrument Sterilization Trays	4 qt	12.75" L X 10.375" W X 2" D
BA054-4	Instrument Sterilization Trays	15 qt	20.75" L X 12.75" W X 4" D
BA054-4H	Instrument Sterilization Trays	7 qt	12.75" L X 10.375" W X 4" D
BA054-6H	Instrument Sterilization Trays	10 qt	12.75" L X 10.375" W X 6" D

MOPEC GROUP

Mopec ELEVATING PATHOLOGY

MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI