SCALES

Digital Compact Scale

Durable, fluid-resistant precision scale for use in all laboratory environments. Offers additional percent weighing mode and security features.



BB063

SPECIFICATIONS

Construction
- LCD Display

Dimensions
- 9.5″ L x 9.8″ W x 2.9″ H

Utilities / Power
- Rechargeable battery, 6 V (100 hour typical operation with 12-hour full charge), or 4 C batteries (not included)

Performance
- Calibration: 1: 20 lb/ 10 kg

Capacity
- 16,000 g

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI