SCALES

Digital Infant Scale

This digital infant scale provides precise weight measurement with a 44 lb capacity and 0.5 oz increments. It includes time/date display and weight-memory functions, with a removable infant tray and AC power adapter.



BB064

SPECIFICATIONS

Construction
- Removable infant tray
- Digital display and control panel

Dimensions
- Scale: 11.3” L X 9.75” W X 7.5” H
- Tray: 10.5” L X 20.6” W X 4” H

Performance
- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- Memory function compares current and stored weight
- Time and date display

Utilities / Power
- (1) 9V battery — sold separately
- (1) 3V CR2032 lithium backup battery
- (1) AC adapter

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI