SCALES

Digital Washdown Platform Scale

This digital washdown platform scale combines durable stainless steel construction with a removable weigh pan for easy cleaning. Its compact design, LED display, and IP68 water- and dust-resistant construction make it well suited for demanding laboratory environments.



BB066

SPECIFICATIONS

Construction
- LED display
- Removable weigh pan

Dimensions
- Platform dimensions: 9.25” X 8”

Performance
- Capacity: 300 oz / 15,000 g
- Readability: 0.1 oz / 5 g
- IP68 water and dust resistant
- Compact and easy to use

Utilities / Power
- AC adapter included

Materials
- Stainless steel

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI