SCALES

Solar Powered Hanging Scale

This solar-powered hanging scale provides accurate digital weighing without the need for replaceable batteries. It includes a hanging pan assembly and supports tare/zero functions for the bow, ring, and pan.



BB067

SPECIFICATIONS

Construction

- Hanging scale
- Includes bow and pan assembly
- Hanging assembly includes:
  Carabiner
  Bow
  Ring and attachment hardware
  Pan

Performance

- Capacity: 30 lb / 15 kg
- Accuracy: 0.02 lb / 0.01 kg
- Digital display provides instant weight readings
- Zero/tare function
- Switchable lb / kg units
- Minimum operating light level: 400 lux

Utilities / Power

- Completely solar powered
- No replacement batteries required

MOPEC GROUP
Mopec
MORTECH MANUFACTURING
www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI